MANUFACTURER Fisher & Paykel Healthcare Ltd






English	А
Español	В
Français	С

TABLE OF CONTENTS

1.	Befo	pre you Start	Α-	1	
2.	Sym	bol Definitions	Α-	1	
3.	Gen	eral Warnings and Cautions	Α-	2	
4.	Inte	nded Use	Α-	3	
5.	Quid	ck Start Guide	Α-	4	
6.	Gett	ing to Know Your F&P ICON™	Α-	5	
	Rep	lacement Parts	Α-	6	
7.	Men	u System	Α-	7	
	(1)	On/Off and Ramp	Α-	7	
	(2)	Humidity	Α-	7	
	(3)	Alarm On/Off	Α-	7	
	(4)	Alarm Time	Α-	7	
	(5)	Alarm Volume	Α-	7	
	(6)	Clock Time	Α-	7	
	(7)	Brightness	Α -	7	
	(8)	Sleep Data	Α -	8	
	(9)	SensAwake™ On/Off	Α -	8	
	(10)	User Preferences	Α -	8	
		(i) Altitude Adjustment	Α -	8	
		(ii) ECO Mode	Α -	8	
		(iii) Display Pressure	Α -	8	
		(iv) Error Code	Α -	9	
8.	Con	nmunicating Sleep Data	Α -	9	
	Tele	phone Reporting	Α -	9	
	Dov	vnloading Data to the SmartStick™	Α-	9	
9.	Clea	aning and Maintenance	Α-	9	
10.	Prod	duct Specifications	Α-	10	
11.	Dev	ice and Consumables Disposal Instructions	Α-	11	
12.	Trou	ubleshooting	Α-	11	
13.	F&P	ICON™ Models and Features Matrix	Α -	11	

1. BEFORE YOU START

Please read this manual carefully before use and save these instructions for future reference.

If you feel that your device or any accessories are not operating correctly, please contact your healthcare provider.

NOTE: This device contains no serviceable parts. Refer all therapy and equipment queries to your healthcare provider.

IMPORTANT CONTACT DETAILS
Home Medical Equipment Provider
Health Insurance Provider

CAUTION: USA Federal law restricts this device for sale by or on the order of a physician.

2. SYMBOL DEFINITIONS

†	Type BF Applied Part	~	Alternating Current
[]i	ATTENTION Consult accompanying documents		Class II Medical Electrical Equipment
X	Do not throw away	€0123	Conforms with medical device directive 93/42/EEC
IPX2	Drip-proof		Date of manufacture
REF	Catalogue number	SN	Serial number
\mathbf{R}_{Only}	Prescription only	EC REP	Authorized representative in the European community
C Z078	NZ Radio Interference C-tick mark	<u>^</u>	Caution

3. GENERAL WARNINGS AND CAUTIONS

This manual refers to the **F&P ICON™** series model as "the device".

The device treats Obstructive Sleep Apnea (OSA) by delivering a flow of continuous positive airway pressure (CPAP) at a level prescribed by the physician, to splint open the airway and prevent airway collapse.

CONTRAINDICATIONS

Research indicates that the following pre-existing conditions may contraindicate
the use of positive pressure for some patients: pneumothorax, bullous lung disease,
pneumocephalus, cerebrospinal fluid leak, recent cranial surgery or trauma,
abnormalities of the cribriform plate, pathologically low blood pressure or in patients
whose upper airways are bypassed.

PRECAUTIONS

- The safety and effectiveness of positive pressure has not been established in patients with respiratory failure or chronic obstructive pulmonary disease.
- The safety and effectiveness of the auto-adjusting device has not been established in patients with congestive heart failure, obesity hypoventilation syndrome or central sleep apnea.

ADVERSE EFFECTS

 Nosebleeds and ear and sinus discomfort may occur from the use of positive pressure therapy.

WARNINGS

To avoid electric shock from the device:

- Only operate if the device, power cord and plug are dry and in good working order.
- If water damage occurs to your device, disconnect the power cord, discontinue use immediately and seek advice from your healthcare provider.
- · Do not store or use the device where it can be pulled into water.

To avoid choking, or inhalation of a foreign object:

- Never place any non-approved objects into any opening of the device, ThermoSmart™ Breathing Tube or mask.
- · Ensure the recommended filter is fitted to the device before use.
- Ensure the ThermoSmart™ Breathing Tube is positioned so it can not become entangled with the body or furniture during sleep.

To avoid burns:

- · Do not fill the chamber with hot water.
- After use, wait for the water to cool before touching, carrying or emptying the Water Chamber.
- · Do not touch the water in the Water Chamber while the device is operating.
- Position the ThermoSmart™ Breathing Tube so it is uncovered and free from bedding or other materials.

To ensure optimal therapy:

- Do not operate the device, Water Chamber or ThermoSmart™ Breathing Tube if it is dropped, damaged or not working as intended.
- Pressure adjustments should only be made by a qualified healthcare provider.
- Only use masks, breathing tubes and accessories compliant with ISO 17510-2, distributed for use with this device, and recommended by Fisher & Paykel Healthcare or your healthcare provider.
- Do not block the exhaust flow holes of the mask as they are designed to allow a continuous flow of air out of the mask and CO₂ re-breathing may occur.
- During a power failure your device will stop running and you will no longer be receiving therapy. If you are using a CPAP mask which meets the specific mask standard (ISO 17510-2), which all Fisher & Paykel Healthcare masks meet, then your CO₂ will be flushed from the mask through the exhaust flow ports. If you wake up and your device has stopped and is unable to be restarted, you should remove your mask and re-initiate therapy once the device has power.
- On models without automatic altitude adjustment, make sure the altitude level is manually adjusted to ensure optimal pressure delivery.
- Do not use the device without a Water Chamber in place.

Other:

- · Refer all repairs and maintenance to your healthcare provider.
- Ensure the device is operated on a level surface, lower than head height, to prevent water entering the tubing and the device enclosure.
- · The device is not intended to be used as a life-supporting device.
- The device complies with the electromagnetic compatibility requirements of IEC 60601-1-2. In certain circumstances the device may affect or be affected by nearby portable mobile radio frequency communication equipment, due to the effects of electromagnetic interference. If this should happen, try moving your device or the location of the equipment causing interference, or alternatively consult your healthcare provider.

Using supplemental oxygen with your device:

- Supplemental oxygen can be administered at the mask end of the breathing tube.
 NOTE: At a fixed flow rate of supplemental oxygen the inhaled oxygen concentration will vary, depending on the pressure settings, the patient's breathing pattern, mask selection and leak rate.
- Ensure there is no obstruction downstream of the Oxygen Port as this can affect the delivered oxygen concentration.
- 3. Oxygen concentration should be measured at the point of delivery to the patient.
- 4. Avoid the risk of fire:
 - Only use oxygen when the device is operating. If the device is turned off it can lead to accumulation of oxygen within the device.
 - · Ensure adequate ventilation is provided around the device.
 - Remove any source of ignition, such as: cigarettes, an open flame, or materials which burn or ignite easily at high oxygen concentration.
 - Keep oxygen regulators, cylinder valves, tubing, connections and all other oxygen
 equipment away from oil, grease or greasy substances. Spontaneous and violent
 ignition may occur if these substances come into contact with oxygen under
 pressure.

CAUTIONS 1

To prevent water damage to the device:

- · Remove the Water Chamber from the device before filling.
- · Do not fill the Water Chamber above the maximum level.
- · Do not move, carry, transport or store the device with water in the Water Chamber.
- · Do not use the device without a Water Chamber.

General:

- Only use the device within the Operating Conditions specified in Section 10.
- Position the device so the power cord connection to the power supply is easily accessible.
- Only clean the device in accordance with the cleaning instructions set out in Section 9 and only when it is disconnected from the power supply.
- Only use the F&P ICON™ SmartStick™ with the device. Use of any other USB drives
 may cause data corruption. Do not attempt to change the directories or view the
 data without software distributed or designed for use with the F&P ICON™.

4. INTENDED USE

The device is for use on adult patients for the treatment of OSA. The device is for use in the home or sleep laboratory.

5. QUICK START GUIDE

Refer to **Part A** of the Diagram Sheet at the beginning of this manual for illustrative examples:

- Place the device on a firm flat surface and plug the power cord into the wall socket.
- Remove the Chamber Lid and reach inside the Water Chamber with your fingers to take hold of the circular chamber grip. Pull the Water Chamber up and out of the device.

CAUTION: Do not fill the Water Chamber while it is in the device; the Water Chamber should always be removed before filling.

- Fill the Water Chamber with distilled water up to the maximum water line, which is marked on the side of the Water Chamber.
 - CAUTION: Do not use the device without a Water Chamber.
- To replace the Water Chamber, the arrow on the circular chamber grip should face toward the Elbow.
- Fit the Chamber Lid by lining up the large dot on the lid with the large dot on the top of the device. When the two symbols are aligned, apply downward force and turn clockwise to align the large dot with the small dot on the top of the device.
- When the Chamber Lid is locked into place, it should be flush, with no gaps for air to leak. It can now be used as a handle for lifting and carrying.
 - NOTE: The Chamber Lid must be correctly locked for optimal pressure delivery.
- 7 The Chamber Lid can be unlocked by twisting in an anti-clockwise direction, then removed by lifting away from the device.
- Attach the breathing tube to the Elbow and connect to your mask.

 NOTE: Make sure the electrical connectors on the ThermoSmart™ Breathing

 Tube come together with the electrical connectors on the Elbow.
- 9 When the selection circle is at the top of the Menu System a quick press to the SmartDial™ will start therapy. The Home Screen will light up and the clock-face of the Menu System will rotate to indicate therapy is starting. If Ramp is required, press the SmartDial™ again and hold for 3 seconds to activate the Ramp. The Ramp symbol will light up on the display.
- Turn the SmartDial™ to select and adjust settings. Refer to Section 7 for a description of the Menu System.

Your device is now ready for use.

6. GETTING TO KNOW YOUR F&P ICON™

Description of Parts and Instructions for Use

Refer to **Part B** at the beginning of this manual for diagrams of each part. Also refer to **Section 5 Quick Start Guide** for additional information.

Water Chamber

 Remove the Water Chamber from the device before filling. Fill to the maximum water line (14.2 fl oz/420 mL), using distilled water only. Do not use hot water.

CAUTION: Do not fill the Water Chamber while it is inside the device and do not fill above the maximum water line.

 To replace the Water Chamber into the device, the arrow on the circular chamber grip should face toward the Elbow.

NOTE: Replace water before each use.

2 SmartDial™

The SmartDial™ is for operating the device and navigating the Menu System. The SmartDial™ works on a simple *Press, Turn, Press* principle. For example:

- Press the SmartDial[™] to turn the device on or off.
- Turn the SmartDial™ to find the required setting.
- Press the SmartDial[™] to enter the setting.
- Turn the SmartDial[™] to adjust the setting.
- Press the SmartDial[™] to accept the change.

NOTE: The SmartDial™ only requires a quick press, unless you are trying to access sub-menus. Sub-menus can be accessed by pressing and holding the SmartDial™ for 3 seconds.

3 The Display

The Display is where information can be viewed. The Home Screen is the default display where each point of the Menu System can be viewed. Once a selection has been made within the Menu System, the Display will dim and return to the Home Screen after 30 seconds.

NOTE: If the device is left idle for 15 minutes the Menu System will lock. To reactivate turn the SmartDial™ 180°.

4 Menu System

The 16 points of the clock-face on the Home Screen also represent each of the available settings.

NOTE: These instructions include functions only available for patient use.

Refer to Section 7 for a detailed description of each setting, and Part C of the Diagram Sheet for illustrative examples.

5 Chamber Lid

- · Twist the Chamber Lid anti-clockwise to unlock and remove.
- The Chamber Lid must be correctly locked into place for optimal pressure delivery.
- Once the Chamber Lid is locked into place, the handle can be used to carry the device.
- The Chamber Lid should be flush with the device with no gaps when it is fitted correctly.

6 ThermoSmart™ Breathing Tube

ThermoSmart™ includes a unique heated breathing tube as well as a humidity algorithm which both work to deliver optimal levels of humidity without the side effects of condensation.

CAUTION: Do not scrunch or crush the ThermoSmart™ Breathing Tube when packing as this can cause damage. If there is any visible damage that exposes copper wires, do not use.

7 Elbow

The Elbow is where the breathing tube connects to the device. The Elbow can be rotated from side to side, depending on your preferred sleeping position.

8 Air Filter

The Air Filter can be found at the rear of the device. To access, first pull out the Filter Grill, then pinch the gauze to remove. The Filter Grill must be pushed back into the cavity with the black line of the filter facing out to keep the Air Filter in place. Replace the Air Filter when it becomes significantly discolored, at least once every three months or after 1,000 hours' worth of running time.

USB Port

- The USB Port houses the SmartStick™ which must be pushed completely into the USB Port until it clicks into place.
- To remove, push the end of the SmartStick™ for a spring-release retrieval.
 NOTE: Only use the SmartStick™ designed and distributed for use with the F&P ICON™.

10 SmartStick™

- The SmartStick™ is used to record sleep data from the device memory when
 it is inserted into the USB Port of the device, and to store music for use with
 AlarmTunes™. Refer to Section 7. Point 3. for more information.
- The SmartStick™ can be removed and sent to your healthcare provider, with the SmartStick™ Mailer to view and adjust settings. Changes will be automatically transferred when the SmartStick™ is next inserted into the USB Port of the device

NOTE: The SmartStick™ must be pulled all of the way out of the USB Port before pushing it back in.

11 Use and Care Manual

Read this manual carefully and save these instructions for future reference. If you have any questions please contact your healthcare provider.

Replacement Parts

900ICON200S Water Chamber (Single) 900ICON214 AUTO Lid

900ICON200M Water Chamber (10 pack) 900ICON216 PREMO/NOVO Lid

900ICON202 SmartStick™ (5 pack) 900ICON218 Filter Grill

900ICON204 Elbow 900ICON315 Soft Carry-case

900ICON206 Outlet Seal

900ICON208 ThermoSmart™ Breathing Tube*

900HC221 Breathing Tube*

900ICON212 Air Filter

* Applied Parts - to fit 0.86 inches (22 mm) Conical Connector

7. MENU SYSTEM

Refer to ${\bf Part}\;{\bf C}$ of the Diagram Sheet at the beginning of this manual for illustrative examples.

1 On/Off and Ramp

- When the selection circle is at the top of the Menu System a quick press to the SmartDial™ will start therapy. The top point of the menu screen will flash to indicate that therapy has commenced.
- To activate Ramp, press the SmartDial™ again and hold for 3 seconds.
 The Ramp symbol will appear .
- To stop and start pressure delivery, press the SmartDial[™] when the selection circle is at the top point.

2 Humidity

To control the level of humidity delivered to the mask turn the SmartDial™ to
the Humidity setting
 I Press to enter and turn the SmartDial™ again to select
the desired humidity level from 0 (off) to 7 (maximum). Press again to accept
the change. For example, if the air is feeling too cold increase the Humidity
setting in increments of 1 or 2 over 2 to 3 nights until you find the level that is
comfortable.

NOTE: The default setting is Humidity level 4.

 The boost control offers three additional levels to optimize the Humidity setting. Press and hold the SmartDial[™] for 3 seconds while in the Humidity setting § and the Boost symbols § § will appear. Turn to select low, medium or high and press again to accept the change.

CAUTION: The Boost control should only be used in extreme situations when experiencing symptoms from your CPAP therapy. The default setting (medium) will be suitable for most patients.

3 Alarm On/Off

- Turn the SmartDial™ to the Alarm setting A, press to enter, then turn to select either On ✓ or Off X. Press again to accept the change.
- If Alarm On

 is chosen, turn the SmartDial™ to select either the default bell

 or the music symbol for AlarmTunes™

 i, and press again to accept the change.
- Press the SmartDial™ once to snooze the Alarm for 10 minutes, twice to turn the Alarm off and discontinue therapy.

NOTE: The AlarmTunes™ symbol is only visible when the SmartStick™ is inserted. SmartStick™ Studio software is required to activate the AlarmTunes™ feature.

 If AlarmTunes™ is selected, turn the SmartDial™ to select the desired music track (from 1 to 5, depending on how many tracks are loaded onto the SmartStick™), or select Auto for shuffle and press again to accept the change. Refer to the SmartStick™ Studio CD for SmartStick™ Studio software, and instructions for installing and transferring music to your SmartStick™.

4 Alarm Time

 Turn the SmartDial™ to the Alarm Time setting ^o⁄₂ and press to enter. Turn the SmartDial™ to select the hour, press to accept. Turn to select the minutes and press again to accept the change.

5 Alarm Volume

- To adjust the volume of the alarm, turn the SmartDial™ to the Alarm Volume setting ⁶⁄₄ and press to enter. Turn the SmartDial™ to adjust the volume from quiet to loud and press again to accept the change.

6 Clock Time

7 Brightness

- Turn the SmartDial™ to the Brightness setting ★ and press to enter. Turn the
 SmartDial™ to select the Home Screen brightness level and press to accept
 the change. To adjust the brightness of the display while accessing the
 Menu System, press and hold the SmartDial™ for 3 seconds and the Active
 Brightness symbol will appear ★ a.
- Turn to select the desired brightness and press again to accept the change.
 NOTE: The Home Screen setting controls the brightness emitted by the LCD when the device is in stand-by mode, and the Active setting controls the brightness emitted by the LCD whilst navigating the Menu System.

8 Sleep Data

Turn the SmartDial^{\mathbb{M}} to the Sleep Data setting $\underline{\mathbb{M}}$ and press to enter. Turn the SmartDial^{\mathbb{M}} to either the CMS report $\underline{\mathfrak{G}}$ (US only) or Patient report $\underline{\mathfrak{ZO}}$ and press to enter.

NOTE: The model name will display while in the Sleep Data setting.

Turn the SmartDial $\!\!\!^{\text{\tiny TM}}$ to scroll through each of the following statistics:

- · Total Compliant Hours.
- Average Compliant Hours over the last 30 nights. Press the SmartDial™ to access the average for the last 7 nights and again for the last night.
- Percentage of nights used greater than 4 hours over the last 30 nights. Press the SmartDial™ to access the average for the last 7 nights.
- The number of nights used since therapy commenced.
- · Check Sum.

NOTE: A ② symbol will appear on the Display Screen when the Sleep Data is within the compliance requirements. Once CMS requirements are met, the data within the CMS sub-menu will freeze.

• Exit the Sleep Data setting () this symbol will only appear if Additional Sleep Data is restricted).

Additional Sleep Data (NOTE: this information may be restricted by your healthcare provider.)

- · AHI for last 30, 7, and last night
- · Leak for last 30, 7, and last night
- 90th percentile pressure for last 30, 7, and last night
- SensAwake[™] for last 30, 7, and last night
- · Exit the Sleep Data setting.

SensAwake™ On/Off

SensAwake™ is activated as a default setting. To deactivate from models with
this feature, turn the SmartDial™ to the SensAwake™ setting ❖ and press to
enter. Turn the SmartDial™ to select On ✔ or Off X and press again to accept
the change.

10 User Preferences 9

- For models without automatic Altitude Adjustment, turn the SmartDial™ to
 the Altitude setting
 <u>M</u> and press to enter. Turn the SmartDial™ to select the
 Altitude level from 0 to 9.000 ft and press to accept the change.
- To choose meters, press the SmartDial™ and hold for 3 seconds while at the Altitude Adjustment setting. When the symbol starts to flash, turn to select meters and press to accept the change. Press the SmartDial™ again to enter the menu, and turn to adjust to the desired Altitude level. Press again to accept the change.

(ii) ECO Mode

ECO Mode allows for the device to operate on a reduced power (75 W).

NOTE: Humidity output will be significantly reduced in this mode and is not recommended for long-term use.

Situations for use in this mode might include long-haul air flights and camping.

- Turn the SmartDial™ to the Eco setting ECO and press to enter. Turn to select
 On ✓ or Off X and press again to accept the change.
- The Eco symbol ECO will be displayed continuously on the Home Screen while the Economy mode is activated.

(iii) Display Pressure

The Display Pressure setting provides two display options for the Home Screen. Time is the default setting, however this can be changed to show pressure (cmH₂O or hPa) if required.

- Turn the SmartDial™ to the Display Pressure setting and press to enter.
 Turn to select On ✓ or Off X, and press again to accept the change.
- If a different unit of measure is required, press the SmartDial™ and hold for 3 seconds. Turn the SmartDial™ to select either cmH₂O or hPa and press again to accept the change.

(iv) Error Code

If the Error symbol *Error* flashes on the Home Screen, turn the SmartDial™ to the User Preferences **(***) and press to enter. Turn the SmartDial™ again to access the Error Code *Error*. Record the number displayed and contact your healthcare provider.

NOTE: Settings will not change when the device is turned off.

8. COMMUNICATING SLEEP DATA

Telephone Reporting

Many healthcare providers will ask for sleep data over the telephone to track CPAP therapy progress. Follow the steps outlined under Sleep Data (Section 7, Point 8) which explains how to access this data.

Downloading Data to the SmartStick™

If your healthcare provider requests a copy of your sleep data from the SmartStick™ please follow the steps below to download this information:

- (1) If the SmartStick™ is already inserted into the USB Port of the device, check that the pressure is off and there is no SmartStick™ symbol ♣ flashing on the Home Screen. The SmartStick™ can now be safely removed by pushing the end for a spring-release retrieval.
- (2) If the SmartStick™ has been removed, ensure the pressure is off before inserting the SmartStick™ back into the USB Port of the device. Wait until the SmartStick™ symbol • ♣ has stopped flashing on the Home Screen before removing the SmartStick™ again.
- (3) Data from the SmartStick™ should then be sent in the manner requested by your healthcare provider.

NOTE: The flashing symbol indicates that data is being transferred from the device memory to the SmartStick™. Removal of the SmartStick™ during this time may cause data corruption.

9. CLEANING AND MAINTENANCE

The device should be cleaned as follows:

- · Unplug the device from the power supply.
- Wipe the exterior and basin of the device (when the Water Chamber has been removed) with a clean, damp (not wet) cloth and mild dishwashing detergent.
 Do not use harsh abrasives, solvents or alcohol as these may damage the device.

Daily

Clean the Water Chamber and the breathing tube:

- · Disconnect the breathing tube from the mask and the device Elbow.
- Clean the breathing tube with warm soapy water, rinse thoroughly and hang with both ends pointing to the floor to dry.
- · Remove the Water Chamber and discard any remaining water.
- The Water Chamber can be washed in a domestic dishwashing machine, or cleaned and washed with mild dishwashing detergent, then rinse and dry.

Weekly

Thoroughly clean the Water Chamber:

- Soak the inside of the Water Chamber for 10 minutes in a solution of one part white vinegar to two parts water.
- Empty the vinegar solution and rinse the Water Chamber thoroughly with water.

It is recommended to replace your Water Chamber and breathing tube regularly. Contact your healthcare provider for advice on how frequently your equipment should be replaced. Equipment should be replaced immediately if there is any sign of cracking, discoloration or leaking.

NOTE: This device is not serviceable and should not be disassembled. Refer all therapy and equipment queries to your healthcare provider.

10. PRODUCT SPECIFICATIONS

DIMENSIONS: 6.3 H x 6.7 W x 8.7 D in. (160 H x 170 W x 220 D mm)

WEIGHT: 4.8 lb (2.2 kg) Packed weight 8.7 lb (4.0 kg)

PERFORMANCE:

Pressure Range: 4 to 20 cmH₂O (in the unlikely event of fault conditions, pressure

may reach up to 26 cmH₂O)

Maximum Flow Rates					
CPAP Pressure Setting (cmH ₂ O)	4	8	12	16	20
Maximum flow rate at patient connection port (L/pM) - AUTO/PREMO	62.1	101.5	129.4	148.8	143.3
Maximum flow rate at patient connection port (L/pM) - NOVO	46.7	52.4	56.5	60.6	64.5

Dynamic Pressure Stability					
CPAP Pressure Setting (cmH ₂ O)	4	8	12	16	20
Pressure difference (cmH ₂ O) at mask connection port - AUTO/PREMO	0.52	0.47	0.47	0.52	0.54
Pressure difference (cmH ₂ O) at mask connection port - NOVO	1.16	1.22	1.28	1.38	1.44

Static Pressure Stability					
	AUTO/PREMO	NOVO			
Change at connection port at a pressure setting of 10 cmH ₂ O	0.2	0.4			

HUMIDITY:

With ThermoSmart™ Breathing Tube:

Maximum humidity = 36 mg/L (BTPS), 82% RH at $10 \text{ cmH}_2\text{O}$, with humidity level 7 and Boost level high. Typical humidity = 24 mg/L (BTPS), 90% RH at $10 \text{ cmH}_2\text{O}$, with humidity level 4 and Boost level medium. With Standard Breathing Tube:

Maximum humidity = 32.24 mg/L, 73.21% RH at 10 cmH₂O, with Humidity Level 7 and Boost Level high. Typical humidity = 17.97 mg/L, 85.88% RH at 10 cmH₂O, with Humidity Level 4 and Boost Level medium.

GAS TEMPERATURES: Maximum = 100 °F (38 °C)

NOISE LEVEL: Sound Pressure Level < 29 dBA
Average Sound Power Level < 37 dBA

WATER CHAMBER VOLUME: 14.2 fl oz (420 mL) up to the maximum fill line

STANDARDS COMPLIANCE: IEC60601-1: 1988 + A1 & A2

AS/NZS 3200.1.0: 1998

EN60601-1: 1990

DATA RECORDING: The SmartStick™ records up to 365 days of summary efficacy data,

30 days of detailed efficacy data (AHI, Leak), Pressure and Leak data (50 Hz), as well as cumulative data from first use.

The device memory can also record up to 365 days of efficacy

data, 400 hours or 200 of detailed efficacy data.

Electrical Ratings		
Rated supply voltage	Rated current input	Rated supply frequency
100 - 115 V	1.27 A (1.43 A Max)	50 - 60 Hz/400 Hz
220 - 240 V	1.07 A (1.21 A Max)	50 - 60 Hz

STORAGE AND TRANSPORT:

The device should always be stored and transported within the

following temperatures:

Transport: 14 to 140 °F (-10 °C to 60 °C) Storage: 14 to 140 °F (-10 °C to 60 °C)

CAUTION: Ensure water is emptied from the Water Chamber

before traveling.

OPERATING CONDITIONS: Ambient Temperature: 41 to 95 °F (5 to 35 °C)

 Humidity:
 5 to 95% Relative Humidity

 Altitude:
 0 to 9,000 ft (0 to 3,000 m)

NOTE: Above 4,500 ft (1,500 m) the maximum operating pressure

will be reduced.

11. DEVICE AND CONSUMABLES DISPOSAL INSTRUCTIONS

Unit Disposal Instructions

This unit contains electronics. Please do not discard as regular rubbish. Dispose of electronics according to local guidelines.

Consumable Disposal Instructions

Place the mask, breathing tube and Water Chamber in a rubbish bag at the end of use and discard with normal rubbish.

12. TROUBLESHOOTING

If you feel that your device is not operating correctly, please refer to the following suggestions. If the problem persists, please consult your healthcare provider. Do not attempt to repair the device yourself.

Problem	Possible Cause	Solution
There is no visible display or power	The plug may not be correctly attached to the power cord.	Ensure the plug is securely attached to the power cord, and is connected to the power supply.
Insufficient air is delivered from the device	The Chamber Lid may not be correctly fitted to the device. Ensure the Water Chamber is in the device.	Refer to Section 5 and Part A for detailed instructions for fitting the Chamber Lid.
"Error" is displayed on the Home Screen of the device	An error may have been detected with the device.	Turn the SmartDial™ to the User Preferences 👣 and press to enter. Turn the SmartDial™ again to access the Error Code Error . Record the number displayed and contact your healthcare provider for further instruction.
AlarmTunes™ does	The Alarm may be set to "off".	Refer to Section 7, Point 3, for instructions on the Alarm
not work	The Alarm may be set to "buzzer".	setting.
SmartStick™ Studio software may not have been installed.		Once SmartStick™ Studio software has been installed onto your computer and your favorite songs selected,
	The Alarm volume may not be set loud enough.	download onto your SmartStick™ and then ensure the SmartStick™ is inserted into the USB Port of the device.
The Heater-Plate and/or Water Chamber base is warm to the touch even though the device isn't being used	The Heater-Plate and Water Chamber are insulated by the device which can cause heat to be retained. This is perfectly normal and should not cause alarm.	If concerned about overheating unplug the device from the wall socket when it is not being used.
I think my chamber is leaking because there is a build-up of water on the Heater-Plate	In stand-by mode the power supply can heat the water in the Water Chamber causing condensation to form.	Before each use remove the Water Chamber and dry the basin of the device with a cloth. If the water build-up becomes excessive, please contact your healthcare provider.
I don't think my humidifier is working	ECO Mode may be activated. When ECO Mode is activated the humidity output will be significantly reduced. Check the Humidity setting is at 4 or above.	To disable ECO Mode turn the dial to the ECO setting ECO and press to enter. Turn to select OFF X and press to accept the change. Refer to Section 7, Point 2, of the manual for more information on the humidity controls.

13. F&P ICON™ MODELS AND FEATURES MATRIX

The model name of each device is located on the base, as well as in the Sleep Data setting of the Menu System. Refer to Section 7, Point 8, for more information.

Features	AUTO	PREMO	NOVO
Fully integrated CPAP with humidification	•	•	
Compatible with ThermoSmart™ Breathing Tube	•	•	•
Auto-Adjusting Pressure	•		
Efficacy Reporting* (AHI, Leak)	•	•	
Compliance Reporting	•	•	•
SmartStick™ Removable Media	•	•	•
SensAwake™ Pressure Relief	•		
Proportional Ramp	•	•	•
Auto-Altitude Adjusting	•	•	Manual
Leak Compensation	•	•	

^{*} NOTE: The display of efficacy data may be restricted by your healthcare provider.

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- Ajuste la tapa de la cámara alineado el punto grande situado en la tapa © con el punto grande situado en la parte superior del dispositivo ©. Cuando se hayan alineado los dos símbolos, presione hacia abajo y gire en el sentido de las agujas del reloj para alinear el punto grande con el pequeño situado en la parte superior del dispositivo.